

— ABLYBIO, Help Your Research



# REST/NRSF Rabbit pAb

货号: **AYP16449**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:100 - 1:500
理论分子量	35kDa/36kDa/119kDa/121kDa
实测分子量	150KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Raji,SH-SY5Y,Mouse heart,Rat heart
细胞定位	Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-130 of human REST/NRSF (NP_005603.3).
------	--

## 靶点信息

研究背景	This gene encodes a transcriptional repressor that represses neuronal genes in non-neuronal tissues. It is a member of the Kruppel-type zinc finger transcription factor family. It represses transcription by binding a DNA sequence element called the neuron-restrictive silencer element. The protein is also found in undifferentiated neuronal progenitor cells and it is thought that this repressor may act as a master negative regulator of neurogenesis. Alternatively spliced transcript variants have been described
基因ID	5978
基因名	REST
Swiss	Q13127 ( <a href="https://www.uniprot.org/uniprotkb/Q13127/entry">https://www.uniprot.org/uniprotkb/Q13127/entry</a> )
别名	REST,NRSF,WT6,XBR,REST/NRSF Rabbit pAb,Neural-restrictive silencer factor,X2 box repressor

## 产品验证



Western blot analysis of REST/NRSF expressed in Mouse heart,Rat heart,HeLa using REST/NRSF Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))