

HOXB7 Rabbit pAb

货号: **AYP16548**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:1000 |
| 理论分子量 | 24kDa |
| 实测分子量 | 24KDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse kidney |
| 细胞定位 | Nucleus |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human HOXB7 (NP_004493.3). |
| 序列 | MSSLYANTLFSKYPASSSVFATGAFPEQTSCAFASNPQRPGYGAGSGASFAASMQGLYPGGGGMAGQSAAGVYAAGY GLEPSSFNMHCAPFEQNLSGVCPGDSAKAAGAKEQRDSDLAA |

靶点信息

| | |
|------|--|
| 研究背景 | This gene is a member of the Antp homeobox family and encodes a protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded nuclear protein functions as a sequence-specific transcription factor that is involved in cell proliferation and differentiation. Increased expression of this gene is associated with some cases of melanoma and ovarian carcinoma. |
| 基因ID | 3217 |

| | |
|-------|---------------------------------|
| 基因名 | HOXB7 |
| Swiss | P09629 |
| 别名 | HOXB7;HHO.C1;HOX2;HOX2C;Hox-2.3 |

产品验证



Western blot analysis of HOXB7 expressed in Mouse kidney using HOXB7 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn