

# DR1 Rabbit pAb

货号: **AYP16559**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	19kDa
实测分子量	19KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis
细胞定位	Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human DR1 (N P_001929.1).
序列	MASSGNDDDLTIPRAAINKMIKETLPNVRVANDARELVVNCCTEFIHLLISSEANEICNKSEKKTISPEHVIAQLES LGFGSYIS EVKEVLQECKTVALKRRKASSRLENLGIPEEELLRQQQLFAKARQQQAELAQQEWLQMQQAAQQAQLAAASASASNQA GSSQDEEDDDI

## 靶点信息

研究背景	This gene encodes a TBP- (TATA box-binding protein) associated phosphoprotein that represses both basal and activated levels of transcription. The encoded protein is phosphorylated in vivo and this phosphorylation affects its interaction with TBP. This protein contains a histone fold motif at the amino terminus, a TBP-binding domain, and a glutamine- and alanine-rich region. The binding of DR1 repressor complexes to TBP-promoter complexes may establish a mechanism in which an altered DNA conformation, together with the formation of higher order complexes, inhibits the assembly of the preinitiation complex and controls the rate of RNA polymerase II transcription.
基因ID	1810
基因名	DR1
Swiss	Q01658
别名	DR1;NC2;NC2-BETA;NC2B;NCB2

## 产品验证



Western blot analysis of DR1 expressed in Mouse testis using DR1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)