

— ABLYBIO, Help Your Research



# IL17F Rabbit pAb

货号: **AYP16708**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	18kDa
实测分子量	18kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,U-87MG
细胞定位	Secreted
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-163 of human IL17F (NP_443104.1).
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a cytokine that shares sequence similarity with IL17. This cytokine is expressed by activated T cells, and has been shown to stimulate the production of several other cytokines, including IL6, IL8, and CSF2/GM-CSF. This cytokine is also found to inhibit the angiogenesis of endothelial cells and induce endothelial cells to produce IL2, TGF $\beta$ 1/TGF $\beta$ , and monocyte chemoattractant protein-1.
基因ID	112744
基因名	IL17F
Swiss	Q96PD4 ( <a href="https://www.uniprot.org/uniprotkb/Q96PD4/entry">https://www.uniprot.org/uniprotkb/Q96PD4/entry</a> )
别名	IL17F,CANDF6,IL-17F,ML-1,ML1,IL17F Rabbit pAb,Cytokine ML-1

## 产品验证



Western blot analysis of IL17F expressed in K-562,U-87MG using IL17F Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)