

TIAL1 Rabbit pAb

货号: **AYP16857**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	41kDa/43kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,K-562
细胞定位	Cytoplasm,Cytoplasmic granule,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 266-375 of human TIAL 1 (NP_003243.1).
序列	EGHVVKCYWGKESPDMTKNFQQVDYSQWGQWSQVYGNPQQYGYMANGWQVPPYGVYQWPWNQQGFGVDQSPS AAWMGGFGAQPPQGQAPPPVIPPNNQAGYGMASYQTQ

靶点信息

研究背景	The protein encoded by this gene is a member of a family of RNA-binding proteins, has three RNA recognition motifs (RRMs), and binds adenine and uridine-rich elements in mRNA and pre-mRNAs of a wide range of genes. It regulates various activities including translational control, splicing and apoptosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. The different isoforms have been shown to function differently with respect to post-transcriptional silencing.
------	---

基因ID	7073
基因名	TIAL1
Swiss	Q01085
别名	TIAL1;TCBP;TIAR

产品验证



Western blot analysis of TIAL1 expressed in A-549, HeLa, K-562 using TIAL1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn