

# SENP2 Rabbit pAb

货号: **AYP16927**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	68kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	rat lung, rat testis
细胞定位	cytosol, nuclear membrane, nuclear pore, nucleoplasm, nucleus, PML body
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 366-589 of human SENP 2 (NP_067640.2).
序列	LELTEDMEKEISNALGHGPQDEILSSAFKLRITRGDIQTLKNYHWLNDEVINFYMNLLVERNKKQGYPALHVFSTFFYPKLKSGGYQAVKRWTKGVNLFQEIILVPIHRKVVHWSLVVIDLRKKCLKYLDSMGQKQHRICEILLQYLQDESKTKRNSDLNLEWTHHSMKPHEIPQQLNGSDCGMFTCKYADYISRDKPITFTQHQMPLFRKKMVWEILHQQLL

## 靶点信息

研究背景	SUMO1 (UBL1; MIM 601912) is a small ubiquitin-like protein that can be covalently conjugated to other proteins. SENP2 is one of a group of enzymes that process newly synthesized SUMO1 into the conjugatable form and catalyze the deconjugation of SUMO1-containing species.[supplied by OMIM, Apr 2004]
基因ID	59343

基因名	SENP2
Swiss	Q9HC62
别名	AXAM2;SMT3IP2;SENP2

## 产品验证



Western blot analysis of SENP2 expressed in rat lung, rat testis using SENP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)