

— ABLYBIO, Help Your Research



IL3 Rabbit pAb

货号: **AYP16969**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	17kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,Mouse thymus
细胞定位	Secreted
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 20-152 of human IL3 (N P_000579.2).
------	---

靶点信息

研究背景	The protein encoded by this gene is a potent growth promoting cytokine. This cytokine is capable of supporting the proliferation of a broad range of hematopoietic cell types. It is involved in a variety of cell activities such as cell growth, differentiation and apoptosis. This cytokine has been shown to also possess neurotrophic activity, and it may be associated with neurologic disorders.
基因ID	3562
基因名	IL3
Swiss	P08700 (https://www.uniprot.org/uniprotkb/P08700/entry)
别名	IL3,IL-3,MCGF,MULTI-CSF,IL3 Rabbit pAb,Hematopoietic growth factor,Mast cell growth factor,Multipotential colony-stimulating factor,P-cell-stimulating factor

产品验证



Western blot analysis of IL3 expressed in Raji, Mouse thymus using IL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)