

— ABLYBIO, Help Your Research



# HEXA Rabbit pAb

货号: **AYP16977**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	19kDa/60kDa
实测分子量	55KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,A-549,HeLa,Mouse liver,Mouse testis,Rat testis,Mouse kidney
细胞定位	Lysosome
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 270-529 of human HEX A (NP_000511.2).
------	---

## 靶点信息

研究背景	This gene encodes a member of the glycosyl hydrolase 20 family of proteins. The encoded preproprotein is proteolytically processed to generate the alpha subunit of the lysosomal enzyme beta-hexosaminidase. This enzyme, together with the cofactor GM2 activator protein, catalyzes the degradation of the ganglioside GM2, and other molecules containing terminal N-acetyl hexosamines. Mutations in this gene lead to an accumulation of GM2 ganglioside in neurons, the underlying cause of neurodegenerative disorders termed the GM2 gangliosidoses, including Tay-Sachs disease (GM2-gangliosidosis type I). Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed.
基因ID	3073
基因名	HEXA
Swiss	P06865 ( <a href="https://www.uniprot.org/uniprotkb/P06865/entry">https://www.uniprot.org/uniprotkb/P06865/entry</a> )
别名	HEXA,TSD,HEXA Rabbit pAb,Beta-N-acetylhexosaminidase subunit alpha,N-acetyl-beta-glucosaminidase subunit alpha

## 产品验证



Western blot analysis of HEXA expressed in Mouse liver,Rat testis,HepG2 using HEXA Rabbit pAb at 1:100. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)