


基因ID	7083
基因名	TK1
Swiss	P04183
别名	TK1;TK2

产品验证



Western blot analysis of Thymidine Kinase 1 expressed in Mouse liver, Rat liver, HeLa using Thymidine Kinase 1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn