

— ABLYBIO, Help Your Research



# Twist Rabbit pAb

货号: **AYP17006**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	<b>IHC:</b> 1:50 - 1:200
理论分子量	20kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Twist (NP_000465.1).
------	---

## 靶点信息

研究背景	Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. The protein encoded by this gene is a bHLH transcription factor and shares similarity with another bHLH transcription factor, Dermo1. The strongest expression of this mRNA is in placental tissue; in adults, mesodermally derived tissues express this mRNA preferentially. Mutations in this gene have been found in patients with Saethre-Chotzen syndrome.
基因ID	7291
基因名	TWIST1
Swiss	Q15672 ( <a href="https://www.uniprot.org/uniprotkb/Q15672/entry">https://www.uniprot.org/uniprotkb/Q15672/entry</a> )
别名	TWIST1,ACS3,BPES2,BPES3,CRS,CRS1,CSO,SCS,TWIST,bHLHa38, Twist Rabbit pAb, Class A basic helix-loop-helix protein 38, H-twist

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)