

— ABLYBIO, Help Your Research



KRT80 Rabbit pAb

货号: **AYP17108**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	47kDa/50kDa/54kDa
实测分子量	47kDa/50kDa/54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3 cells, human skeletal muscle tissue
细胞定位	cytoplasm,cytosol,intermediate filament cytoskeleton
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 242-452 of human KRT80 (NP_872313.2).
------	---

靶点信息

研究背景	Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into epithelial keratins and hair keratins. This gene's expression profile shows that it encodes a type II epithelial keratin, although structurally the encoded protein is more like a type II hair keratin. This protein is involved in cell differentiation, localizing near desmosomal plaques in earlier stages of differentiation but then dispersing throughout the cytoplasm in terminally differentiating cells. The type II keratins are clustered in a region of chromosome 12q13. Two transcript variants encoding two different fully functional isoforms have been found for this gene.
基因ID	144501
基因名	KRT80
Swiss	Q6KB66 (https://www.uniprot.org/uniprotkb/Q6KB66/entry)
别名	KRT80,KB20,KRT80 Rabbit pAb,Cytokeratin-80,Keratin-80,Type-II keratin Kb20

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)