

# NKRF Rabbit pAb

货号: **AYP17134**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	77kDa/79kDa
实测分子量	78KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,Mouse brain,Mouse spleen,Mouse kidney,Rat brain
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 391-690 of human NKRF (NP_060014.2).
序列	YGTKKTSKHAADAALKILQKTQPTYPSVKSSQCHTGSSPRGSGKKKDIKDLVYENSSNPVCTLNDAQFNRMTVEYVY ERMTGLRWKCKVILESEVIAEAVGVKKTVKYEAAGEAVKTLKKTQPTVINNLKKGAVEDVISRNEIQGRSAEEAYKQKIKED NIGNQLLRKMGWTGGGLGKSGEGIREPISVKEQHKREGLDVERVKNIAKRDIEQIIRNYARSESHTDLTFSRELTNDERK QIHQIAQKYGLKSKSHGVGHDRYLVVGRKRKEDLLDQLKQEGQVGHYELVMPQAN

## 靶点信息

研究背景	This gene encodes a transcriptional repressor that interacts with specific negative regulatory elements to mediate transcriptional repression of certain nuclear factor kappa B responsive genes. The protein localizes predominantly to the nucleolus with a small fraction found in the nucleoplasm and cytoplasm. Alternative splicing results in multiple transcript variants.
------	--

基因ID	55922
基因名	NKRF
Swiss	O15226
别名	NKRF;ITBA4;NRF

## 产品验证



Western blot analysis of NKRF expressed in Mouse brain,Rat brain,U-87MG using NKRF Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)