

— ABLYBIO, Help Your Research



# UXT Rabbit pAb

货号: **AYP17220**

## 产品信息

|       |  |
|-------|--|
| 反应    | Human  |
| 宿主    | Rabbit   |
| 克隆性   | Polyclonal   |
| 预测反应  |  |
| 应用    | WB   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000  |
| 理论分子量 | 18kDa  |
| 实测分子量 | 17kDa  |
| 形式    | Liquid   |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated   |
| 阳性对照  | 293T,HL-60   |
| 细胞定位  | Cytoplasm,Nucleus,centrosome,cytoskeleton,microtubule organizing center                                |
| 纯化    | Affinity purification  |


## 抗原信息

|      |  |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 13-169 of human UXT (NP_705582.1). |
|------|--|

## 靶点信息

|       |  |
|-------|--|
| 研究背景  | The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. |
| 基因ID  | 8409   |
| 基因名   | UXT  |
| Swiss | Q9UBK9 ( <a href="https://www.uniprot.org/uniprotkb/Q9UBK9/entry">https://www.uniprot.org/uniprotkb/Q9UBK9/entry</a> )   |
| 别名    | UXT,ART-27,STAP1,UXT Rabbit pAb,Androgen receptor trapped clone 27 protein,Ubiquitously expressed transcript protein   |

## 产品验证



Western blot analysis of UXT expressed in 293T,HL-60 using UXT Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)