

— ABLYBIO, Help Your Research



EXTL3 Rabbit pAb

货号: **AYP17266**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:200 - 1:2000 IHC: 1:50 - 1:200
理论分子量	104kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse pancreas,Rat lung
细胞定位	Endoplasmic reticulum membrane,Single-pass type II membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 52-220 of human EXTL3 (NP_001431.1).
------	--

靶点信息

研究背景	This gene encodes a single-pass membrane protein which functions as a glycosyltransferase. The encoded protein catalyzes the transfer of N-acetylglucosamine to glycosaminoglycan chains. This reaction is important in heparin and heparan sulfate synthesis. Alternative splicing results in the multiple transcript variants.
基因ID	2137
基因名	EXTL3
Swiss	O43909 (https://www.uniprot.org/uniprotkb/O43909/entry)
别名	EXTL3,BOTV,EXTL1L,EXTR1,REGR,RPR,ISDNA,EXTL3 Rabbit pAb,EXT-related protein 1,Glucuronyl-galactosyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase,Hereditary multiple exostoses gene isolog,Multiple exostosis-like protein 3,Putative tumor suppressor protein EXTL3

产品验证



Western blot analysis of EXTL3 expressed in 293T, Mouse pancreas, Rat lung using EXTL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)