

THBS3 Rabbit pAb

货号: **AYP17293**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	91kDa/104kDa
实测分子量	104kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse ovary
细胞定位	extracellular region,perinuclear region of cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 787-956 of human THBS 3 (NP_009043.1).
序列	TDDDYAGFLFSYQDSGRFYVVMWKQTEQTYWQATPFRAVAQPGLQLKAVTSVSGPGEHLRNALWHTGHTPDQVRLW TDPNVGWRDTSYRWQLLRPQVGYIRVKLYEGPQLVADSGVIIDTSMRGGRLGVFCFSQENIIWSNLQYRCNDTVPE DFEPFRQLLQGRV

靶点信息

研究背景	The protein encoded by this gene belongs to the thrombospondin family. Thrombospondin family members are adhesive glycoproteins that mediate cell-to-cell and cell-to-matrix interactions. This protein forms a pentameric molecule linked by a single disulfide bond. This gene shares a common promoter with metaxin 1. Alternate splicing results in coding and non-coding transcript variants.
------	--

基因ID	7059
基因名	THBS3
Swiss	P49746
别名	THBS3;TSP3

产品验证



Western blot analysis of THBS3 expressed in Mouse kidney, Mouse ovary using THBS3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn