

# U2SURP Rabbit pAb

货号: **AYP17318**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:20 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	72kDa/118kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,LO2,Jurkat
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 705-1029 of human U2SURP (NP_001073884.1).
序列	GAPLEDVDGIPIDATPIDDLDGVPIKSLDDDDLDGVPLDATEDSKKNEPIFKVAPSKWEAVDESELEAQAVTTSKWELFDQH EESEEEENQNQEESEDEEDTQSSKSEEHHLYSNPIKEEMTESKFSKYSEMSEEKRAKLREIELKVMKFQDELESGKRPKPK GQSFQEQVEHYRDKLLQREKEKELERERERDKKDKKESRSKDKKEDECTPTRKERKRRHSTSPSPSRSSSGRRVKSP SPKSERSERSERSHKESSRSRSSHKDSPRDVSKKAKRSPSGSRTPKRSRRRSRSPKKSRSRSRSPHRSRSHKKSCKN KH

## 靶点信息

研究背景	Enables RNA binding activity. Predicted to be involved in RNA processing. Located in nucleoplasm.
------	---

基因ID	23350
基因名	U2SURP
Swiss	O15042
别名	U2SURP;SR140;fSAPa

## 产品验证



Western blot analysis of U2SURP expressed in HeLa, LO2, Jurkat using U2SURP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)