

— ABLYBIO, Help Your Research



# UBE4A Rabbit pAb

货号: **AYP17323**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:1000 - 1:5000 <b>IHC:</b> 1:50 - 1:200
理论分子量	122kDa/123kDa
实测分子量	123KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Cytoplasm
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 797-1066 of human UBE 4A (NP_001191006.1).
------	--

## 靶点信息

研究背景	This gene encodes a member of the U-box ubiquitin ligase family. The encoded protein is involved in mult i ubiquitin chain assembly and plays a critical role in chromosome condensation and separation through t he polyubiquitination of securin. Autoantibodies against the encoded protein may be markers for sclerode rma and Crohn's disease. A pseudogene of this gene is located on the long arm of chromosome 3. Alterna tively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	9354
基因名	UBE4A
Swiss	Q14139 ( <a href="https://www.uniprot.org/uniprotkb/Q14139/entry">https://www.uniprot.org/uniprotkb/Q14139/entry</a> )
别名	UBE4A,E4,UBOX2,UFD2,UBE4A Rabbit pAb,RING-type E3 ubiquitin transferase E4 A,KIAA0126

## 产品验证



Western blot analysis of UBE4A expressed in Mouse liver,Rat liver using UBE4A Rabbit pAb at 1:1000. Sec ondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buf fer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)