

UBE4A Rabbit pAb

货号: **AYP17323**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:1000 - 1:5000 IHC: 1:50 - 1:200
理论分子量	122kDa/123kDa
实测分子量	123KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 797-1066 of human UBE 4A (NP_001191006.1).
序列	FLLDEAIQYLSKIKIQIEKDRGEWDSLTPPEARREKEAGLQMFQQLARFHNIMSNETIGTLAFLTSEIKSLFVHPFLAERIISMLNYFLQHVLVGPKMGALKVKDFSEDFKPKQLVSDICTIYLNLGDEENFCATVPKDGSRYSPTLFAQTVRVLKKINKPGNMIMAFSNLAERIKSLADLQQQEEETYADACDFLDPIMSTLMCDPVVLPSSRVTVDRSTIARHLLSDQTDPFNRSPLTMDQIRPNTLKEKIQRWLAERKQKQEQLE

靶点信息

研究背景	This gene encodes a member of the U-box ubiquitin ligase family. The encoded protein is involved in mult ubiquitin chain assembly and plays a critical role in chromosome condensation and separation through the polyubiquitination of securin. Autoantibodies against the encoded protein may be markers for scleroderma and Crohn's disease. A pseudogene of this gene is located on the long arm of chromosome 3. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	--

基因ID	9354
基因名	UBE4A
Swiss	Q14139
别名	UBE4A;E4;UBOX2;UFD2

产品验证



Western blot analysis of UBE4A expressed in Mouse liver, Rat liver using UBE4A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn