




基因ID	7936
基因名	NELFE
Swiss	P18615
别名	NELFE;D6S45;NELF-E;RD;RDBP;RDP

### 产品验证



Western blot analysis of NELFE expressed in Mouse liver, HeLa, 293T using NELFE Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)