

MCMBP Rabbit pAb

货号: **AYP17391**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	52kDa/72kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,Mouse spleen,Mouse liver,Mouse ovary
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human MCMBP (NP_079110.1).
序列	MPCGEDWLSHPLGIVQGFFAQNGVNPDWEKKVIEYFKEKLNENAPKWVPSLNEVPLHYLKPNSFVKFRCMIQDMFDPE FYMGVYETVNQNTKAHVLHFGKYRDVAECGPQQELDLNSPRNTTLERQTFYCVVPGESTWVKEAYVNNANQARVSPSTS YTPSRHKRSYEDDDMDLQPNKQKQDHAGARQAGSVGGLQWCGEPKRLETEASTGQQLNSLNLSSPFDLNFPLPGEK GPACL

靶点信息

研究背景	This gene encodes a protein which is a component of the hexameric minichromosome maintenance (MCM) complex which regulates initiation and elongation of DNA. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	79892

基因名	MCMBP
Swiss	Q9BTE3
别名	MCMBP;C10orf119;MCM-BP

产品验证



Western blot analysis of MCMBP expressed in Mouse spleen, A-431, Mouse liver using MCMBP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn