

# FAF1 Rabbit pAb

货号: **AYP17423**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	56kDa/73kDa
实测分子量	85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR3,HT-29,Mouse brain,Mouse liver,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human FAF1 (N P_008982.1).
序列	MASNMDREMILADFQACTGIENIDEAITLLEQNNWDLVAAINGVIPQENGILQSEYGGGETIPGPAFNPASHPASAPTSSSSS AFRPVMPSRQIVERQPRMLDFRVEYRDRNVDVLEDTCTVGEIKQILENELQIPVSKMLLKGWKTGDVEDSTVLKSLHLPK NNSLYVLTPDLPPSSSSSHAGALQESLNQNFMLIITHREVQREYNLNFSGSSTIQEVKRNVDLTSIPVRHQLWEGWPTSA TDDSMCLAESGLSYPCHRLTVGRRSSPAQTREQSEE

## 靶点信息

研究背景	Interaction of Fas ligand (TNFSF6) with the FAS antigen (TNFRSF6) mediates programmed cell death, also called apoptosis, in a number of organ systems. The protein encoded by this gene binds to FAS antigen and can initiate apoptosis or enhance apoptosis initiated through FAS antigen. Initiation of apoptosis by the protein encoded by this gene requires a ubiquitin-like domain but not the FAS-binding domain.
------	--

基因ID	11124
基因名	FAF1
Swiss	Q9UNN5
别名	FAF1;CGI-03;HFAF1s;UBXD12;UBXN3A;hFAF1

## 产品验证



Western blot analysis of FAF1 expressed in Mouse brain,Rat brain,OVCAR3 using FAF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)