

— ABLYBIO, Help Your Research



SIAH1 Rabbit pAb

货号: **AYP17522**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC IF/ICC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 |
| 理论分子量 | 21kDa/31kDa/34kDa |
| 实测分子量 | 31kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse kidney |
| 细胞定位 | Cytoplasm,Nucleus |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 140-282 of human SIAH 1 (NP_003022.3). |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | This gene encodes a protein that is a member of the seven in absentia homolog (SIAH) family. The protein is an E3 ligase and is involved in ubiquitination and proteasome-mediated degradation of specific proteins. The activity of this ubiquitin ligase has been implicated in the development of certain forms of Parkinson's disease, the regulation of the cellular response to hypoxia and induction of apoptosis. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized. |
| 基因ID | 6477 |
| 基因名 | SIAH1 |
| Swiss | Q8IUQ4 (https://www.uniprot.org/uniprotkb/Q8IUQ4/entry) |
| 别名 | SIAH1,SIAH1A,SIAH1 Rabbit pAb,RING-type E3 ubiquitin transferase SIAH1,Seven in absentia homolog 1,Siah-1a,HUMSIAH |

产品验证



Western blot analysis of SIAH1 expressed in Mouse kidney using SIAH1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)