

ARL3 Rabbit pAb

货号: **AYP17585**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	20kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,LO2,HeLa,Mouse eye,Mouse brain,Mouse heart,Mouse kidney,Rat brain
细胞定位	Cell projection,Cytoplasm,Cytoplasmic side,Golgi apparatus membrane,Nucleus,Peripheral membrane protein,centrosome,cilium,cytoskeleton,microtubule organizing center,spindle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of human ARL3 (NP_004302.1).
序列	MGLLSILRKLKSAPDQEVRIILLGLDNAGKTTLLKQLASEDISHITPTQGFNIKSVQSQGFKLNVDIGGQRKIRPYWKNYF ENTDILYVIDSADRKRFEETGQELAEELLEEEKLSCVPVLIFANKQDLLTAAPASEIAEGLNLHTIRDRVWQIQSCSALTGEGV QDGMNVVCKNVNAKKK

靶点信息

研究背景	ADP-ribosylation factor-like 3 is a member of the ADP-ribosylation factor family of GTP-binding proteins. A RL3 binds guanine nucleotides but lacks ADP-ribosylation factor activity.
------	---

基因ID	403
基因名	ARL3
Swiss	P36405
别名	ARL3;ARFL3

产品验证



Western blot analysis of ARL3 expressed in Mouse eye,Rat brain,A-431 using ARL3 Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn