

— ABLYBIO, Help Your Research



SAA4 Rabbit pAb

货号: **AYP17630**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	14kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,Mouse spleen,Rat heart
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human SAA4 (NP_006503.2).
------	--

靶点信息

研究背景	Major acute phase reactant.
------	-----------------------------

基因ID	6291
基因名	SAA4
Swiss	P35542 (https://www.uniprot.org/uniprotkb/P35542/entry)
别名	SAA4,C-SAA,CSAA,SAA4 Rabbit pAb,Constitutively expressed serum amyloid A protein

产品验证



Western blot analysis of SAA4 expressed in Mouse spleen,Rat heart,293T using SAA4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)