

# PCF11 Rabbit pAb

货号: **AYP17699**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	173kDa
实测分子量	180kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Jurkat,HeLa,Mouse testis
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1256-1555 of human PC F11 (NP_056969.2).
序列	PKAYPDNHLNQVDVNELEFSKLLKLGILKLSQTDSATTQVSEVTAQPPPEEEEDQNEQDQVDPDLTNTVEELKQRYDSVINR LYTGIQCYSCGMRFTTSQTDVYADHLDWHYRQNRTEKDVSRKVTHRRWYYSLTDWIEFEEIADLEERAKSQFFEKVHEEV VLKTQEAAKEKEFQSVAPAGPAGAVESCEICQEQFEQYWDEEEEEEWHLKNAIRVDGKIYHPSCYEDYQNTSSFDCTPSPSK TPVENPLNIMLNIVKNELQEPDCKVKEERIDTPPACTEESIATPSEIKTENDTVESV

## 靶点信息

研究背景	The protein encoded by this gene binds to CLP1 to form pre-mRNA cleavage factor II <sub>m</sub> . The encoded protein is necessary for efficient Pol II transcription termination and may be involved in degradation of the 3' product of polyA site cleavage.
基因ID	51585

基因名	PCF11
Swiss	O94913
别名	PCF11

## 产品验证



Western blot analysis of PCF11 expressed in Mouse testis, LO2, Jurkat using PCF11 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)