

# ZCCHC10 Rabbit pAb

货号: **AYP17704**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	18kDa/20kDa
实测分子量	21kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye
细胞定位	
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human ZCCHC 10 (NP_001287745.1).
序列	MATPMHRLIARRQAFDTELQPVKTFWILIQPSIVISEANKQHVRCQKCLEFGHWTYECTGKRKYLHRPSRTAELKKALKEKENRLLLQSIGETNVERKA

## 靶点信息

研究背景	Predicted to enable nucleic acid binding activity and zinc ion binding activity.
基因ID	54819
基因名	ZCCHC10

Swiss	Q8TBK6
别名	ZCCHC10

## 产品验证



Western blot analysis of ZCCHC10 expressed in Mouse eye using ZCCHC10 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)