

— ABLYBIO, Help Your Research



# ZNF471 Rabbit pAb

货号: **AYP17730**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	28kDa/73kDa
实测分子量	73KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG, HEL, Molt-4 (Negative control)
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 80-210 of human ZNF471 (NP_065864.2).
------	---

## 靶点信息

研究背景	May be involved in transcriptional regulation.
------	--

基因ID	57573
基因名	ZNF471
Swiss	Q9BX82 ( <a href="https://www.uniprot.org/uniprotkb/Q9BX82/entry">https://www.uniprot.org/uniprotkb/Q9BX82/entry</a> )
别名	ZNF471,ERP1,Z1971,ZNF471 Rabbit pAb,EZFIT-related protein 1,KIAA1396

## 产品验证



Western blot analysis of ZNF471 expressed in U-251MG,HEL,Molt-4(Negative control) using ZNF471 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)