

# TNKS2 Rabbit pAb

货号: **AYP17733**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	126kDa
实测分子量	126kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse liver,mouse kidney,rat liver
细胞定位	Chromosome,Cytoplasm,Golgi apparatus membrane,Nucleus,Peripheral membrane protein,telomere
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 780-1050 of human TNK S2 (NP_079511.1).
序列	TPLDLVSADDVSALLTAAMPPSALPSCYKQVLNGVRSFGATADALSSGPSSPSSLSAASSLDNLSGFSFSELSVSSSGT EGASSLEKKEVPGVDFSTQFVRNLGLEHLMDFEREQITLDVLVEMGHKELKEIGINAYGHRHKLKIGVERLISGQQGLNPY LTLNTSGSGTILIDLSPDDKEFQSVEEEMQSTVREHRDGGHAGGIFNRYNILKIQKVCNKKLWERYTHRRKEVSEENHNHA NERMLFHGSPFVNAIHKGFDERHAY

## 靶点信息

研究背景	Poly-ADP-ribosyltransferase involved in various processes such as Wnt signaling pathway, telomere length and vesicle trafficking. Acts as an activator of the Wnt signaling pathway by mediating poly-ADP-ribosylation of AXIN1 and AXIN2, 2 key components of the beta-catenin destruction complex: poly-ADP-ribosylated target proteins are recognized by RNF146, which mediates their ubiquitination and subsequent degradation. Also mediates poly-ADP-ribosylation of BLZF1 and CASC3, followed by recruitment of RNF146 and subsequent ubiquitination. Mediates poly-ADP-ribosylation of TERF1, thereby contributing to the regulation of telomere length. Stimulates 26S proteasome activity.
基因ID	80351
基因名	TNKS2
Swiss	Q9H2K2
别名	TNKS2;ARTD6;PARP-5b;PARP-5c;PARP5B;PARP5C;TANK2;TNKL;pART6

## 产品验证



Western blot analysis of TNKS2 expressed in mouse liver, mouse kidney, rat liver using TNKS2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)