

— ABLYBIO, Help Your Research



ZGPAT Rabbit pAb

货号: **AYP17739**

产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	20kDa/54kDa/55kDa/57kDa
实测分子量	57KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 212-511 of human ZGPA T (NP_852150.2).
------	--

靶点信息

研究背景	Transcription repressor that specifically binds the 5'-GGAG[GA]A[GA]A-3' consensus sequence. Represses transcription by recruiting the chromatin multiprotein complex NuRD to target promoters. Negatively regulates expression of EGFR, a gene involved in cell proliferation, survival and migration. Its ability to repress genes of the EGFR pathway suggest it may act as a tumor suppressor. Able to suppress breast carcinoma.
基因ID	84619
基因名	ZGPAT
Swiss	Q8N5A5 (https://www.uniprot.org/uniprotkb/Q8N5A5/entry)
别名	ZGPAT,GPATC6,GPATCH6,KIAA1847,ZC3H9,ZC3HDC9,ZIP,ZGPAT Rabbit pAb,G patch domain-containing protein 6,Zinc finger CCCH domain-containing protein 9,Zinc finger and G patch domain-containing protein

产品验证



Western blot analysis of ZGPAT expressed in Mouse spleen using ZGPAT Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)