

DHX57 Rabbit pAb

货号: **AYP17744**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	49kDa/73kDa/155kDa
实测分子量	150-155kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Rat eye
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 700-900 of human DHX 57 (NP_945314.1).
序列	SATLNAELFSDYFNSCPVITIPGRTFPVDQFFLEDAIAVTRYVLQDGSPYMRSMKQISKEKLKARRNRRTAFEEVEEDLRLSLH LQDQDSVKDAVPDQQLDFKQLLARYKGVSKSVIKTMSIMDFEKVNLELIEALLEWIVDGKHSYPPGAILVFLPGLAEIKMLY EQLQNSLSLFNRRRSNRCVHPLHSSLSSEEQAVFV

靶点信息

研究背景	Probable ATP-binding RNA helicase.
基因ID	90957
基因名	DHX57

Swiss	Q6P158
别名	DHX57;DDX57

产品验证



Western blot analysis of DHX57 expressed in U-87MG,Rat eye using DHX57 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer : 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn