

— ABLYBIO, Help Your Research



ZNF763 Rabbit pAb

货号: **AYP17761**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	31kDa/46kDa
实测分子量	46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 315-370 of human ZNF763 (NP_001012771.1).
------	-------------------------------------------------------------------------------------------------------------------------

靶点信息

研究背景	May be involved in transcriptional regulation.
------	------------------------------------------------

基因ID	284390
基因名	ZNF763
Swiss	Q0D2J5 (https://www.uniprot.org/uniprotkb/Q0D2J5/entry)
别名	ZNF763,ZNF,ZNF440L,ZNF763 Rabbit pAb

产品验证



Western blot analysis of ZNF763 expressed in HepG2, Mouse liver using ZNF763 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)