

FCN2 Rabbit pAb

货号: **AYP17802**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:100 - 1:200 |
| 理论分子量 | 34kDa |
| 实测分子量 | 34kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | mosue spleen,mouse lung,rat spleen |
| 细胞定位 | blood microparticle,external side of plasma membrane,extracellular exosome,extracellular region,extracellular space |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 230-290 of human FCN2 (NP_004099.2). |
| 序列 | SAGDSLTFHNNQSFSTKDQDNDLNTGNCAVMFQGAWWYKNCHVSNLNGRYLRGTHGSFANG |

靶点信息

| | |
|------|---|
| 研究背景 | The product of this gene belongs to the ficolin family of proteins. This family is characterized by the presence of a leader peptide, a short N-terminal segment, followed by a collagen-like region, and a C-terminal fibrinogen-like domain. This gene is predominantly expressed in the liver, and has been shown to have carbohydrate binding and opsonic activities. Alternatively spliced transcript variants encoding different isoforms have been identified. |
|------|---|

| | |
|-------|--------------------------------|
| 基因ID | 2220 |
| 基因名 | FCN2 |
| Swiss | Q15485 |
| 别名 | FCN2;EBP-37;FCNL;P35;ficolin-2 |

产品验证



Western blot analysis of FCN2 expressed in mouse spleen, mouse lung, rat spleen using FCN2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn