

— ABLYBIO, Help Your Research



# OTUD6A Rabbit pAb

货号: **AYP17822**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	34kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant protein of human OTUD6A
------	-------------------------------------

## 靶点信息

研究背景	Deubiquitinating enzyme that hydrolyzes 'Lys-27'-, 'Lys-29'- and 'Lys-33'-linked polyubiquitin chains. Also able to hydrolyze 'Lys-11'-linked ubiquitin chains (By similarity).
------	---

基因ID	408193
基因名	Otud6a
Swiss	Q6IE21 ( <a href="https://www.uniprot.org/uniprotkb/Q6IE21/entry">https://www.uniprot.org/uniprotkb/Q6IE21/entry</a> )
别名	OTUD6A,DUBA2,HSHIN6,EG408193,OTUD6A Rabbit pAb,Hin-6 protease,Hin6

## 产品验证



Western blot analysis of OTUD6A expressed in Mouse testis,Rat testis using OTUD6A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)