

ARSA Rabbit pAb

货号: **AYP17988**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	44kDa/53kDa
实测分子量	66kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,HepG2,A-549,Mouse testis,Mouse liver,Mouse brain
细胞定位	Lysosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 210-509 of human ARSA (NP_000478.3).
序列	LMADAQRQDRPFLLYASHHHTYPQFSGQSFAERSGRGPFGDSLMEldaAVGTLMTAIGDLGlleETLVIFTADNGPETM RMSRGGCSGLLRcGKGTTYEGGVREPALAFWPGHIAPGVTHELASSLDLLPTLAALAGAPLPNVTLDGFDLSPLLLGTGKS PRQSLFFYPSYPDEVRGVFAVRTGKYKAHFFTQGSahSDTTADPACHASSSLTAHEPPLLYDLSKDPGENYNLLGGVAGA TPEVLQALKQLLKAQLDAAVTFGSPQVARGEDPALQICCHPGCTPRPACCHCPDPHA

靶点信息

研究背景	The protein encoded by this gene hydrolyzes cerebroside sulfate to cerebroside and sulfate. Defects in this gene lead to metachromatic leucodystrophy (MLD), a progressive demyelination disease which results in a variety of neurological symptoms and ultimately death. Alternatively spliced transcript variants have been described for this gene.
------	---

基因ID	410
基因名	ARSA
Swiss	P15289
别名	ARSA;MLD

产品验证



Western blot analysis of ARSA expressed in Mouse testis, Jurkat, HeLa using ARSA Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn