


基因ID	5077
基因名	PAX3
Swiss	P23760
别名	PAX3;CDHS;HUP2;WS1;WS3

产品验证



Western blot analysis of PAX3 expressed in 293T,Rat testis using PAX3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn