

ZNF620 Rabbit pAb

货号: **AYP18027**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa/48kDa
实测分子量	35kDa/48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human ZNF620 (NP_001243096.1).
序列	MVGGLPGNVSQHLDFGSSLEQPQGHWIIKTKSKRRHFTDTSARHHEAYEVKNGEKFEKLGKNISVSTQLTTNQTNPSTGQI SYECGQCGRYFIQMADFHRHEKCHTGEKSFCEKCEGKYFRYNLLIRHQIIHTGKKPFKCKEKGKGLSSD

靶点信息

研究背景	May be involved in transcriptional regulation.
基因ID	253639
基因名	ZNF620

Swiss	Q6ZNG0
别名	ZNF620

产品验证



Western blot analysis of ZNF620 expressed in mouse lung using ZNF620 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn