

— ABLYBIO, Help Your Research



MARCH6 Rabbit pAb

货号: **AYP18043**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	90kDa/97kDa/102kDa
实测分子量	102kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Endoplasmic reticulum membrane,Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 10-80 of human MARCH 6 (NP_005876.2).
------	---

靶点信息

研究背景	This gene encodes a member of a family of membrane-associated E3 ubiquitin ligases containing RING-C H-type zinc finger motifs. Ubiquitination of type II deiodinase by the encoded protein is an important regulatory step in thyroid hormone signalling. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	10299
基因名	MARCH6
Swiss	O60337 (https://www.uniprot.org/uniprotkb/O60337/entry)
别名	MARCH6,DOA10,MARCH-VI,RNF176,TEB4,MARCH6 Rabbit pAb,Doa10 homolog,Membrane-associated RING finger protein 6,Membrane-associated RING-CH protein VI,Protein TEB-4,RING finger protein 176,RING-type E3 ubiquitin transferase MARCHF6,KIAA0597,44261

产品验证



Western blot analysis of MARCH6 expressed in Mouse brain,Rat brain using MARCH6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)