

ACTL9 Rabbit pAb

货号: **AYP18082**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	45kDa
实测分子量	46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-416 of human ACTL9 (NP_848620.3).
序列	MDASRPKSSSEQSSLEAPRPGPNPSPNVVNKPLQRDSPGMVADRLPPKTGAVVIDMGTGTCKVGFAGQASPTYTVATILGCQPKKPATSGQSLQTFIGEAAARVLPCLTLVQPLRSGIVVDWDAELIWRHLLHDLRVATHDHPFLFSDPPFPATNREKLVEVAFESLRSPAMYVASQSVLSVYAHGRVSGLVVDTGHGVITYVPVFQGYNLLHATERLDLAGNNLTAFLAEMLLQAGLPLGQQDLDLVENIKHHYCYVASDFQKEQARPEQYKRTLKLPDGRVTTLGKELFQCELLFNPPEVPGLSPVGLSTMAKQSLRKLKLSLEMRADLAQNVLLCGSSLFTGFEGRFRAELLRALPAETHVVVAAQPTRNFSVWIGGSILASLRAFQSCWVLRQYEEQGPYIVYRKCY

靶点信息

研究背景	Predicted to be located in cytoplasm. Predicted to be part of dynactin complex. Implicated in spermatogenic failure 53.
------	---

基因ID	284382
基因名	ACTL9
Swiss	Q8TC94
别名	ACTL9;HSD21

产品验证



Western blot analysis of ACTL9 expressed in Mouse testis,Rat testis using ACTL9 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn