

TIGD1 Rabbit pAb

货号: **AYP18085**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	67kDa
实测分子量	67kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse liver,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TIGD1 (NP_663748.1).
序列	MASKCSSERKSRTSLTLNQLKLEMIKLSEEGMSKAEIGRRLGLLRQTVSQVVNAKEKFLKEVKSATPMNTRMIRKRNSLIAD MEKVLVVWIEDQTSRNIPLSQSLIQNKALTLFNSMKAERGVEAAEEKFEASRGWFMRFKERSHFHNIKAQGEAASADVEA AASYPEALAKIIDEGGYTKQIFNVDETAFYWKKMPSRTFIAREEKSVPGFKASKDRLTLLLGANAAGDFKLPMLIYHSEN PRALKNYTKSTLPVLYKWNSKARMTAHLFTAWFTEYFKPTVETYCSEKKIPFKILL

靶点信息

研究背景	The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known.
------	---

基因ID	200765
基因名	TIGD1
Swiss	Q96MW7
别名	TIGD1;EEYORE

产品验证



Western blot analysis of TIGD1 expressed in Mouse brain, Mouse liver, Rat liver using TIGD1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn