

# ARMC6 Rabbit pAb

货号: **AYP18114**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	51kDa/54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,293T,LO2,Mouse brain,Mouse liver,Mouse pancreas,Rat testis,Rat brain
细胞定位	cytosol
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ARMC6 (NP_001186125.1).
序列	MSERCCSRYSSGASIGCTPTSTQAKMVSKRIAQETFDAAVRENIEEFAMGPPEEAVKEAVEQFESQGVDSLNIKTAPKVSADGSQEPHTDILQMLSDLQESVASSRPQEVSAYLTRFCDQCKQDKACRFLAAQKGAYPIIFTAWKLATAGDQGLLLQSLNALSVLTDGQPDLDDAQLQLLVATLTQNADEADLTCSGIRCVRHACLKHEQNRQDLVKAG

## 靶点信息

研究背景	The function of this gene's protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	93436

基因名	ARMC6
Swiss	Q6NXE6
别名	ARMC6;R30923_1

## 产品验证



Western blot analysis of ARMC6 expressed in Mouse brain,Rat testis,A-431 using ARMC6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)