

# USP42 Rabbit pAb

货号: **AYP18131**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	144kDa/145kDa
实测分子量	170kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Raji,SGC7901,Mouse testis
细胞定位	cytosol,nucleoplasm,nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human USP4 2 (NP_115548.1).
序列	ASVQNWSVNRSSVIPEHPKKQKISISHNKLPVRQCQSQPNLHSNSLENPTKPVPSSTITNSAVQSTSNASTMSVSSKVTKP IPRSESCSQPVMNGKSKLNSSVLVPYGAESSEDSDEESKGLGKENGIGTIVSSHSPGQDAEEDDEATPHELQEPMTLNGAN SADSDDPKENGLAPDGASCQGGQPALHSENPFAKANGLP

## 靶点信息

研究背景	Deubiquitinating enzyme which may play an important role during spermatogenesis.
基因ID	84132
基因名	USP42

Swiss	Q9H9J4
别名	USP42

## 产品验证



Western blot analysis of USP42 expressed in Mouse testis, HT-29, Raji using USP42 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)