

FAM98A Rabbit pAb

货号: **AYP18212**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	55kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,LO2,A-549
细胞定位	cytoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human FAM98 A (NP_056290.3).
序列	MECDLMETDILESLEDLGYKGPLLEDGALSQAVSAGASSEFTKLCAWLVSELRVLCKLEENVQATNSPSEAEFQLEVSG LLGEMNCPYLSLTSGDVTKRLLIQKNCLLLTYLISELEAARMLCVNAPPKKAQEGGGSEVFQELKGICIALGMSKPPANITM FQFFSGIEKKLKETLAKVPPNHVGKPLLKPPMPGAHWEKIEAINQAIANEYEVRRK

靶点信息

研究背景	Positively stimulates PRMT1-induced protein arginine methylation. Involved in skeletal homeostasis (By similarity). Positively regulates lysosome peripheral distribution and ruffled border formation in osteoclasts (By similarity). Promotes colorectal cancer cell malignancy.
基因ID	25940

基因名	FAM98A
Swiss	Q8NCA5
别名	FAM98A

产品验证



Western blot analysis of FAM98A expressed in U-87MG, LO2, A-549 using FAM98A Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn