

# REXO2 Rabbit pAb

货号: **AYP18213**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:100
理论分子量	23kDa/26kDa
实测分子量	27kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-87MG,LO2,HT-29
细胞定位	Cytoplasm,Mitochondrion intermembrane space,Mitochondrion matrix,Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 26-237 of human REXO 2 (NP_056338.2).
序列	VREGGAAMAAGESMAQRMVWVDLEMTGLDIEKDQIEMACLITDSDLNILAEGPNLIKQPDELLDSMSDWCKEHHGKSG LTKAVKESTITLQQAEYEFVSFVRQQTTPGLCPLAGNSVHEDKKFLDKYMPQFMKHLHYRIIDVSTVKELCRRWYPEEYFA PKKAASHRALDDISESIKELQFYRNNIFKKKIDKRRKIENGNEKTVS

## 靶点信息

研究背景	This gene encodes a 3'-to-5' exonuclease specific for small (primarily 5 nucleotides or less in length) single-stranded RNA and DNA oligomers. This protein may have a role in DNA repair, replication, and recombination, and in RNA processing and degradation. It may also be involved in resistance of human cells to UV-C-induced cell death through its role in the DNA repair process.
------	---

基因ID	25996
基因名	REXO2
Swiss	Q9Y3B8
别名	REXO2;CGI-114;REX2;RFN;SFN

## 产品验证



Western blot analysis of REXO2 expressed in 293T,U-87MG,LO2 using REXO2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)