

Icam5 Rabbit pAb

货号: **AYP18267**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	96kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 586-786 of mouse Icam 5 (NP_032345.2).
序列	VEGSGKLFSCVEVDGKPEPRVECVGSEGASEGIVLPLVSSNSGPRNSMTPGNLSPGIYLCNATNRHGSTVKTVVVSASPP QMDESSCPHQWLEGAEATALACSARGRPSRVHCSREGAARLERLQVSREDAGTYRCVATNAHGDSRTVTVGVEY RPVVAELAASPPSVRPGGNFTLTCRAEAWPPAQISWRAPPAL

靶点信息

研究背景	ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2).
基因ID	15898
基因名	

Swiss	Q60625
别名	ICAM5;TLCN;TLN;Icam5

产品验证



Western blot analysis of Icam5 expressed in Mouse brain using Icam5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn