


Swiss	P43268
别名	ETV4;E1A-F;E1AF;PEA3;PEAS3;ETS variant 4;Pea3/ETV4

产品验证



Western blot analysis of ETV4 expressed in K-562,Rat brain using ETV4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn