

Septin 12 Rabbit pAb

货号: **AYP18335**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:200 - 1:2000 |
| 理论分子量 | 35kDa/40kDa |
| 实测分子量 | 41kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse brain,Rat brain |
| 细胞定位 | Cell projection,Cytoplasm,cilium,cytoskeleton,flagellum,spindle |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 189-358 of human Septin 12 (NP_653206.2). |
| 序列 | VVPVIARADSLTMEEREAFRRRIQQNLRTHCIDVYPQMCFDEDINDKILNSKLRDRIPFAVVGADQEHLVNGRCVLGRKTK WGIIEVENMAHCEFPLLRDLLIRSHLQDLKDITHNIHYENYRVIRLNESHLLPRGPGWVNLAPASPGQLTTPRTFKVCRGAH DSDDEF |

靶点信息

| | |
|------|--|
| 研究背景 | This gene encodes a guanine-nucleotide binding protein and member of the septin family of cytoskeletal GTPases. Septins play important roles in cytokinesis, exocytosis, embryonic development, and membrane dynamics. Multiple transcript variants encoding different isoforms have been found for this gene. |
| 基因ID | 124404 |

| | |
|-------|-------------------------|
| 基因名 | SEPT12 |
| Swiss | Q8IYM1 |
| 别名 | SEPT12;SPGF10;septin-12 |

产品验证



Western blot analysis of Septin 12 expressed in Mouse brain,Rat brain using Septin 12 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn