

MAF1 Rabbit pAb

货号: **AYP18359**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	28kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,Mouse thymus
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-256 of human MAF1 (NP_115648.2).
序列	MKLENSSEAINSQLTVETGDAHIIGRIESYSCKMAGDDKHMFKQFCQEGQPHVLEALSPPQTSGLSRSLKSKSQGGEEE GPLSDKCSRKTLFYLIATLINESFRPDYDFSTARSHSREPSLSWVNAVNCSLFSAVREDFKDLKPQLWNAVDEEICLAE CDIYSYNDLSDPFGEDGSLWSFNFFYNKRLKRIVFFSCRSISGSTYTPSEAGNELDMELGEEVEEESRSGGSGAEET STMEEDRVPVICI

靶点信息

研究背景	This gene encodes a protein that is similar to Maf1, a <i>Saccharomyces cerevisiae</i> protein highly conserved in eukaryotic cells. Yeast Maf1 is a negative effector of RNA polymerase III (Pol III). It responds to changes in the cellular environment and represses pol III transcription. Biochemical studies identified the initiation factor TFIIIB as a target for Maf1-dependent repression.
------	--

基因ID	84232
基因名	MAF1
Swiss	Q9H063
别名	MAF1;MAF1 homolog

产品验证



Western blot analysis of MAF1 expressed in A-549, HeLa, Mouse thymus using MAF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn