

PCNP Rabbit pAb

货号: **AYP18400**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	6kDa/13kDa/18kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human PCNP (NP_065090.1).
序列	MADGKAGDEKPEKSQRAGAAGGPPEEEAEKPVKTKTVSSSNGGESSRSAEKRSAEAAAADLPTKPTKISKFGFAIGSQT KKASAISIKLGSKPKETVPTLAPKTLVAAAFNEDEDESEPEEMPEAKMRMKNIGRDTPTSAGPNSFNKKGKHGFS DNQKL WERNIKSHLGNVHDQDN

靶点信息

研究背景	May be involved in cell cycle regulation.
基因ID	57092
基因名	PCNP

Swiss	Q8WW12
别名	PCNP

产品验证



Western blot analysis of PCNP expressed in Mouse liver,Rat liver using PCNP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn