

— ABLYBIO, Help Your Research



# ZNF331 Rabbit pAb

货号: **AYP18412**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:200 - 1:2000
理论分子量	53kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-130 of human ZNF331 (NP_061025.5).
------	---

## 靶点信息

研究背景	This gene encodes a zinc finger protein containing a KRAB (Kruppel-associated box) domain found in transcriptional repressors. This gene may be methylated and silenced in cancer cells. This gene is located within a differentially methylated region (DMR) and shows allele-specific expression in placenta. Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding the same protein.
基因ID	55422
基因名	ZNF331
Swiss	Q9NQX6 ( <a href="https://www.uniprot.org/uniprotkb/Q9NQX6/entry">https://www.uniprot.org/uniprotkb/Q9NQX6/entry</a> )
别名	ZNF331,RITA,ZNF361,ZNF463,ZNF331 Rabbit pAb,C2H2-like zinc finger protein rearranged in thyroid adenomas,Zinc finger protein 361,Zinc finger protein 463

## 产品验证



Western blot analysis of ZNF331 expressed in Mouse testis,Rat testis using ZNF331 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))