

— ABLYBIO, Help Your Research



# UTP18 Rabbit pAb

货号: **AYP18432**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:200 - 1:2000
理论分子量	62kDa
实测分子量	60-70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,Mouse brain,Mouse kidney,Mouse liver
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 70-260 of human UTP18 (NP_057085.2).
------	--

## 靶点信息

研究背景	Involved in nucleolar processing of pre-18S ribosomal RNA.
------	--

基因ID	51096
基因名	UTP18
Swiss	Q9Y5J1 ( <a href="https://www.uniprot.org/uniprotkb/Q9Y5J1/entry">https://www.uniprot.org/uniprotkb/Q9Y5J1/entry</a> )
别名	UTP18,CGI-48,WDR50,UTP18 Rabbit pAb,WD repeat-containing protein 50

## 产品验证



Western blot analysis of UTP18 expressed in Mouse brain, HeLa, Jurkat using UTP18 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)